


Oxystat 100/200 Upgrade

GCS 670



Gas Control Systems, Inc. now has the long awaited update for the Oxystat 100 and 200. This GCS 670 offers all the popular features available in the GCS 610 and the GCS 660 that have been available for the past two years.

The 670 uses a PLC that takes the place of the datagem. Pictured left, it mounts on the door of an Oxystat 200 or in a wall box in the case of the Oxystat 100. The PLC operator interface is fully functional and can be programmed directly or through a Windows PC frontend.

Added features the update offers are:

- 640 Temperature Interface updates temperatures immediately, not just each hour with CA readings.
- Flash Memory Card retains data if PC is disconnected for a few days.
- PC logs hourly temps and CA readings. Easily Graphed or Printed
- Temperature Probe allocation and averaging available.
- Temperature Control easily included.
- Multiple systems can communicate down one data line. GCS 610, GCS 660, and now GCS 670 can be linked to a single PC.
- Mexican Protocol accepted with added features and formatted printouts for export record keeping.
- Easily Installed.



420 S. State St. Sparta, MI 49345
ph. 616 887-7994 fax 616 887-1128
800-487-7994
www.gascontrolsystems.com

Distributor: